
Reichert Model 13104190/13104230

Reichert Typ 13104190/13104230

*Automatically Temperature Compensated
Hand-Held Refractometers*

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**Automatically
Temperature Compensated
Hand-Held Refractometers, Reichert
Models 13104190 and 13104230**

Instruction Manual

1.0 Introduction

Reichert Model 13104190 and 13104230 Automatically Temperature Compensated Hand-Held Refractometers are specifically designed for rapid, accurate and reliable readings of total dissolved solids in aqueous solutions. These refractometers can be used to take readings by hand or dipped in a sample solution in a laboratory, factory or field. Temperature compensation occurs automatically.

These instructions will guide you through the proper use and care of Reichert Model 13104190 and 13104230 Automatically Temperature Compensated Hand-Held Refractometers. They also describe techniques and troubleshooting methods to maximize the efficiency of the instruments.

2.0 Measuring Scales

The Reichert Model 13104190 Automatically Temperature Compensated Hand-Held Refractometer measures refractive index (n_D). Full-scale range is 1.3330 to 1.3730. For salinity, full-scale range is 0-160 parts per thousand.

The Reichert Model 13104230 Automatically Temperature Compensated Hand-Held Refractometer measures degrees Brix (% solids) on a full-scale range from 0° to 30° Brix. The smallest scale division is 0.1. All optics and instrument scales are permanently sealed inside the refractometer body to prevent damage.

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3.0 Operating Instructions

Both Reichert Models 13104190 and 13104230 refractometers provide accuracy and reproducibility comparable to bench-type Abbe refractometers. They are automatically temperature compensated for temperatures ranging from 60 ° to 100 F (16 ° to 38 C). The prism is permanently sealed to enable use of the refractometers as dipping instruments.

Follow these steps for immediate, direct accurate readings.

1. Hold the refractometer in a horizontal position (Figure 1). Lift the cover plate (A) to expose both the prism and cover plate surfaces. Clean the cover plate using distilled water. For best results always make sure the prism is clean before placing a sample on it.

Using a plastic pipette or medicine dropper, place a small drop of the sample on the measuring prism (B). Avoid touching the pipette tip to the prism surface.

Immediately close the cover plate over the measuring prism to minimize evaporation.

An alternate method of sample preparation (Figure 2) uses capillary action. After cleaning the prism surface with distilled water and drying it with a soft cloth, close the cover plate. Next, place a drop of the sample liquid on the exposed area on top or bottom of the prism area. The liquid will be drawn into the space between the prism by capillary action. Take care to avoid lifting the cover plate before the reading is taken.

2. To take a direct reading, make sure that the instrument's cover plate is in direct contact with the prism surface. Point the refractometer toward a window or other light source such as a lamp. Aim the refractometer toward the light source until you can see a distinct contrast between the light and dark boundary on the instrument scale (Figures 3 and 4).

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
Figure 3 Model 10419
Abb. 3 Typ 10419


Figure 4 Model 10423
Abb. 4 Typ 10423

3. Focus the scale by rotating the eyepiece. As long as the same individual continues to use the instrument, the focus does not have to be readjusted.

4. Take the reading where the dividing line crosses between the light and dark areas.

5. To clean the prism, use a soft cloth or lens paper moistened with distilled water. Wipe it dry with a soft cloth.

If the prism surface or cover plate are not cleaned completely before the sample is loaded, an erroneous or fuzzy reading will result.

Do not immerse the eyepiece or black focusing ring in water. Never use hot water or abrasive cleaning compounds on the external coating or prism surface. They can damage the instrument.

NOTE: Never expose the instrument to temperatures above 150° F (66° C).

4.0 Calibration

For calibration purposes, your instrument needs to be kept at a temperature between 70° and 85° F (21° and 29° C). To assure proper instrument calibration, take a reading with distilled water.

If the reading departs from zero by more than .05%, gently pry through the cement prism seal and turn the adjustment screw slightly clockwise to increase reading (Fig. 5) with the appropriate tool. Make sure that the final motion is clockwise. Take another reading.

Repeat this procedure until the adjustment is properly made to obtain a valid reading.

Reseal the hole after a correct reading has been obtained.


Figure 5 / Abb. 5

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5.0 Temperature Compensation

Reichert Model 13104190 and 13104230 hand-held refractometers are automatically temperature compensated to provide accurate, direct readings of aqueous solutions at temperatures from 60° to 100° F (16° to 38°C). Maximum temperature compensation error for temperature ranges and for the instrument extremes is 0.1%. The actual error over the most frequently useful portion of the range is much less (see Figure 6).

The largest error occurs at the lowest scale readings and particularly with water samples at temperatures below 65°F (18° C). If water is used for calibration, the temperature should be between 70° and 80°F (21° - 29°C).

Reichert Model 13104190 and 13104230 hand-held refractometers provide automatic temperature compensation required for aqueous solutions. Solid standards are measured accurately only at temperatures between 65° and 70°F (18° and 21°C).

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6.0 Dipping Refractometry

Both Reichert Models 13104190 and 13104230 hand-held refractometers can be used to take readings as dipping instruments. Certain pre-cautions, however, must be taken to assure accuracy.

When used as hand-held refractometers, samples will take on the temperature of the instrument almost instantaneously because the sample is small. However, if the refractometer is used as a dipping instrument it should be allowed to rest in the sample solution three to five minutes, especially if the temperature of the sample differs substantially from that of the refractometer.

The refractometer prism must be in temperature equilibrium with the sample in order to obtain a sharp shadow line. If the temperature equilibrium is poor, a sharp shadowline will not be produced.

7.0 Air Bubble

An air bubble, trapped at the end of the refractometer cavity, accommodates the thermal expansion of the liquid sample.

During transit or under severe vibration, the bubble may escape the trap and appear in the visible portion of the refractometer prism. If this occurs, hold the refractometer vertically, eyepiece facing down, and shake it lightly. This will allow the bubble to pass back into the trap where it will remain during normal operation.

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8.0 Non-Aqueous Solutions

At temperatures between 65° and 75°F (18° and 21°C), readings will generally be correct for all products. The temperature coefficient of refractive index (nD) of oils, hydrocarbons or other liquid products is generally larger than that of water.

Precautions should be taken if materials other than aqueous solutions are to be measured with Reichert Models 13104190 and 13104230 hand-held refractometers.

9.0 Calibrating for Non-Aqueous Solutions

The Brix scale (% solids) is applicable for readings of aqueous solutions provided that the specific product is calibrated to the refractometer.

To calibrate, select three to five accurately known samples (e.g., 2-4-6-8-10%) that have been carefully measured and mixed, preferably by weight (e.g., 2-4-6-8-10%). Determine their Brix (% solids) values using the hand-held refractometer.

Plot a chart with refractometer readings on the vertical y-axis and known Brix (% solids) values on the horizontal x-axis. Unknown solutions can be read in % concentration to an accuracy of 0.1%.

Because Reichert Models 13104190 and 13104230 hand-held refractometers are temperature compensated, correction for this variable is unnecessary.

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©2003 Reichert, Inc. Depew NY, USA 10423-101 Rev. F
AMETEK GmbH