
# Reichert Brix Scale Refractometer

# Reichert Brix Skala Refraktometers

Automatically Temperature Compensated  
Hand-held Refractometers  
Models 137530L0/137531L0/137533L0

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**Automatically Temperature  
Compensated Handheld  
Refractometers  
Models 137530L0/137531L0/137533L0**

**INSTRUCTION MANUAL**

**1.0 INTRODUCTION**

Temperature compensated handheld refractometers from Reichert (Models 137530L0/137531L0/137533L0) are specifically designed for rapid, accurate and reliable readings of total dissolved solids in aqueous solutions. Temperature compensation occurs automatically for samples from 0° to 104°F (-18° to 40°C). Refractometers can easily be used in factory, laboratory or field environments.

These instructions describe proper use and are of the instrument. Techniques and troubleshooting methods to maximize the efficiency

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

**Figure 1 Model 137530L0**  
**Abb. 1 Typ 137530L0**

**Figure 2 Model 137531L0**  
**Abb. 2 Typ 137531L0**

of the instrument are also described.

## 2.0 MEASURING SCALES

The Brix scale originated in the food industry where "Brix" is the primary unit of measurement for dissolved solids. It corresponds to the % sucrose concentration in a solution (e.g., 5°Brix = 5% wt./wt. sucrose).

The refractive index to concentration relationship differs from product to product. This makes it necessary to calibrate the refractometer to specific products and establish a calibration chart for reference (See 7.0 - Calibrating for Other Aqueous Solutions).

**Model 137530L0** - % Solids (°Brix) readings 0° to 30° range. The scale reads directly in % solids (°Brix) with the smallest increment being 0.2° (Figure 1).

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readings 0° to 50° range. The scale reads directly in % Solids (°Brix) with the smallest increment being 0.25° (Figure 2).

**Model 137533L0** - % Solids (°Brix) readings 50° to 90° range. The scale reads directly in % Solids (°Brix) with the smallest increment being 0.2° (Figure 3).

### 3.0 OPERATING INSTRUCTIONS

Automatically temperature compensated hand-held refractometers provide accurate, reproducible readings of total dissolved solids in aqueous solutions.

They are automatically temperature compensated for temperatures ranging from 0° to 104°F (-18° to 40°C). The measuring prism is permanently sealed to enable use of the refractometers as dipping instruments.

Follow these steps for immediate, direct readings:

1. Lift the sample cover plate to expose the measuring prism. Clean the sample cover and measuring prism using distilled water.

Using the dipstick or dropper applicator, place a small drop of the sample on the measuring prism. Avoid touching the pipette tip to the prism surface. Close the sample cover plate

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Figure 4 / Abb. 4


Figure 5 / Abb. 5

immediately to minimize evaporation.

An alternate method of sample application is shown in Figure 4. After cleaning the measuring prism surface with distilled water and drying it with a soft cloth, close the sample cover. Next, place a drop of the sample liquid on the exposed area at the top of the measuring prism area. The liquid will be drawn into the space between the measuring prism and the sample cover plate. Avoid lifting the sample cover plate before the reading is taken.

2. To take a direct reading, the sample cover must be in direct contact with the measuring prism surface. Point the refractometer toward a window or other light source, such as a lamp (Figure 5). A distinct line (where the field changes from light to dark) can be seen on the instrument scale (Figures 1 - 3). Take the reading at this line.

3. After each use thoroughly clean the measuring prism with distilled water and wipe it dry with a soft cloth. An erroneous or fuzzy reading will result if the instrument is not thoroughly cleaned before each use.

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Figure 6 / Abb. 6

water. Hot water or abrasive cleaning compounds can severely damage the external coating or measuring prism surface of the refractometer. Never expose the instrument to temperatures above 150°F (66°C).

#### 4.0 CALIBRATION

Reichert's automatically temperature compensated hand-held refractometers are factory calibrated.

To assure proper readings, the refractometer should be kept at a temperature between 70° and 85°F (21° and 29°C).

To check calibration of models **137530L0** and **137531L0**, take a reading with distilled water. If the reading departs from zero by more than one scale division, pry through the label covering the calibration access hole on the underside of the instrument to expose the phillips head calibration screw (Figure 6). Turn the adjustment screw slightly clockwise to increase the reading. Take another reading.

Repeat this procedure until the adjustment allows a valid reading to be obtained.

After the reading is taken, seal the hole made by the screwdriver point with a sealant such as silicone caulk.

**Model 137533L0** cannot be calibrated with water due to the fact that the instrument is designed for solutions having a minimum concentration of 50% sucrose (1.4201 refractive index equivalent).

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To calibrate this instrument, prepare a solution of a known concentration by weight of sucrose vs water (minimum of 50%). Thoroughly clean the prism with water, wipe dry with a clean, soft cloth. Read the solution on the instrument. If the reading differs from the known value by more than 0.2°Brix, use a jeweler's screwdriver to pry through the seal covering the calibration access hole on the label side of the instrument. Adjust the calibration screw until the reading matches the known value. When complete, reseal the hole with a silicone/RTV sealant to maintain watertight performance.

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
Figure 7 / Abb. 7


Figure 8 / Abb. 8

If assistance is required, please contact our customer service department in the United States at 1 716-686-4500.

### **5.0 TEMPERATURE COMPENSATION**

Automatically temperature compensated hand-held refractometers provide accurate, direct readings of aqueous solutions at temperatures from 0° to 104°F (-18° to 40°C). The temperature of the sample has little bearing on the accuracy of the reading. The sample is so small that it immediately assumes the temperature of the instrument (Figure 7).

### **6.0 NON-AQUEOUS SOLUTIONS**

At temperatures between 65° and 75° (18° and 24°C) readings will generally be correct for all products. The temperature coefficient of refractive index (ne) of oils, hydrocarbons or other liquid products is generally larger than that of water. These samples should be read as close to 68°F (20°C) as possible.

### **7.0 CALIBRATING FOR OTHER AQUEOUS SOLUTIONS**

The % Solids (°Brix) scale is applicable for readings of aqueous solutions provided that the refractometer is calibrated to the specific product.

To calibrate, select three to five accurately known samples (e.g., 2-4-6-8-10%) that have been carefully measured and mixed, preferably by weight determine their % Solids (°Brix) values by using the hand-held refractometer (Figure 8).

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Plot a chart with refractometer readings on the vertical y-axis and known % Solids (°Brix) values on the horizontal x-axis. Unknown solutions can be read in % concentration to an accuracy of 0.1% (Figure 8).

Because these hand-held refractometers automatically temperature compensate, correction for this variable is unnecessary.

**Reichert Analytical Instruments  
Reichert, Inc.**

3362 Walden Ave.  
Depew, NY 14043  
USA

Toll Free: 888-849-8955  
Phone: 716-686-4500  
Fax: 716-686-4555  
email: [ai@reichert.com](mailto:ai@reichert.com)  
[www.reichert.com](http://www.reichert.com)

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Hubertusstrasse 4  
D-82229 Seefeld  
Germany  
Tel: +49-8152-993530

ISO-9001 Certified

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**Reichert GmbH**